


EXPERIMENT 1

AIM: To verify Lambert-beer's law for KMnO_4 colorimetrically.

THEORY: The Beer-Lambert law states that the absorbance of a solution is directly proportional to the concentration of the absorbing species in the solution and the path length. Thus, for a fixed path length (cuvette length), UV/Vis spectroscopy can be used to determine the concentration of the absorber in a solution. The absorbance changes with concentration, A higher concentration of the colored solution absorbs more light (and transmits less) than a solution of lower concentration.

According to Beer-Lambert law,



$$\log (I_o / I_t) = A = \epsilon c l$$

where

I_o and I_t are the incident and transmitted intensities,

A = **absorbance** and

ϵ is a constant *i.e.* **absorptivity** (formerly called the *extinction coefficient*). If the concentration is measured in molL^{-1} , the absorptivity is called the *molar absorptivity*.

$$A = \epsilon c l$$

At constant length

$$A \propto c$$

REQUIREMENTS

Colorimeter, cuvette, six test tubes

Two Burettes or graduated cylinders two 100 mL beakers

0.01M KMnO_4 solution

distilled water, test tube rack, stirring rod, tissues (preferably lint-free)

PROCEDURE

(a) Determination of λ_{\max}

(b) Absorbance of different concentration solution at λ_{\max}

(a) Calculation of λ_{\max}

This can be taken from references (tables of molar extinction coefficients), or more accurately, determined from a calibration curve.

- i. Prepare 100 ml, 0.01M **KMnO₄** (**Molecular weight 158.03 gm/mol**) solution (stock solution), fill it in a burette
- ii. Switch on the computer **and/or** the instrument powers; wait for 30 minutes for ‘warm-up’ of the instrument.
- iii. In the instrument one can choose % transmittance or absorbance display, wavelength range of interest, etc.
- iv. Take one clean and dry glass (only for visible range scan) or quartz cuvette with a given path length (say, 1 cm). Prepare a *blank* by filling a cuvette 3/4 full with distilled water

Note: To correctly use cuvettes, remember:

- Wipe the outside of each cuvette with a lint-free tissue.
 - Handle cuvettes only by the top edge of the ribbed sides.
 - Dislodge any bubbles by gently tapping the cuvette on a hard surface.
 - Always position the cuvette so the light passes through the clear sides.
- v. By using blank cuvette (with distilled water) calibrate the colorimeter (absorbance = 0, transmittance = 100 %) at **filter 1** and put out the blank cuvette.
 - vi. Then fill another cuvette with the stock solution and measure absorbance, note the reading in table.
 - vii. Repeat the step v and vi for all the filters (**wave length range**)
 - viii. From table we can see the λ_{\max} *i.e.* wave length at which solution show maximum absorbance/ O.D.


Observation Table for calculating the λ_{\max}

S.No	Filter	Wavelength/ Range	O.D./ Absorbance
1	1	420	
2	2	440	
3	3	490	
4	4	520	
5	5	540	
6	6	570	
7	7	600	
8	8	650	
9	9	700	

(b) Absorbance of different concentration solutions

- ✓ Obtain small volumes 0.01M KMnO_4 (**Molecular weight 158.03 gm/mol**) of solution and distilled water in separate beakers, fill in the separate graduated burettes
- ✓ Label five clean, dry, test tubes 1–5. Use Burettes to prepare five standard solutions according to the chart below. Thoroughly mix each solution with a stirring rod. Clean and dry the stirring rod between uses.

Concentration can be calculate by $M_1V_1 = M_2V_2$




Test Tube	0.01M KMnO ₄ . (mL)	Distilled H ₂ O (mL)	Concentration(M)
1	2	8	0.002
2	4	6	0.004
3	6	4	0.006
4	8	2	0.008
5	~10	0	0.0100

- ✓ You are now ready to collect absorbance-concentration data for the five standard solutions.
- ✓ Switch on the computer **and/or** the instrument powers; wait for 30 minutes for 'warm-up' of the instrument.
- ✓ In the instrument one can select light sources (UV and visible), choose the slit width, scan speed and % transmittance or absorbance display, wavelength range of interest, etc.
- ✓ Take two clean and dry glass (only for visible range scan) or quartz cuvettes with a given path length (say, 1 cm). Prepare a *blank* by filling a cuvette 3/4 full with distilled water and the other cuvette with aqueous KMnO₄ solution with lowest concentration.
- ✓ Read the absorbance value displayed in the meter. When the displayed absorbance value stabilizes, record its value in your data table.
- ✓ Repeat the procedure for Test Tubes 2 to 5. Similarly spectral runs are done for all the other samples starting from the lowest concentrations to next higher concentrations of KMnO₄. Every time one should rinse the cuvette taking a small portion of the solution to be analyzed next.
- ✓ Plot a curve between Absorbance v/s concentrations. Check whether it is a liner plot or not

Data Table/ Observation Table


Test Tube	0.01M KMnO ₄ . (mL)	Distilled H ₂ O (mL)	Concentration(M)	Absorbance
	0	10	0.00	
1	2	8	0.002	
2	4	6	0.004	
3	6	4	0.006	
4	8	2	0.008	
5	~10	0	0.0100	



Result A linear curve is obtained between Absorbance v/s concentrations that prove the existence of Lambert-Beer law.

PRECAUTIONS


One should note that the Beer–Lambert law is obeyed by many substances mainly at low to moderate concentrations; therefore, dilute concentrations of the absorbing species should be measured. In practice it is advisable to measure absorbances in the range $0.1 < A < 1.0$. Care must be taken to avoid any kind of chemical associations/dissociations of the absorbing species.




Colorimeter



UV-Visible Spectrophotometers



Cuvette



Instrumentation

SPECTROPHOTOMETER

