MOHANLAL SUKHADIA UNIVERSITY, UDAIPUR B. Sc. BIOTECHNOLOGY III YEAR TDC (2016-17)

B. Sc. IIIrd Year Practicals

Practical - I

(A) Microbial Technology

- 1. Isolation of industrially important microorganisms for microbial processes.
- 2. To test the production of enzymes: Amylase, proteinases, lipases and celluloses by microorganisms.
- 3. Demonstration of citric acid production by *Aspergillus niger*, *Penicillium citrianum*.
- 4. To study industrial production of beer/ wine.
- 5. Demonstration of production of antibiotics (penicillin) by microbes.
- 6. To study general methods of food preservation (e.g. Temperature, Salt, Moisture).
- 7. Testing of milk by MBRT.
- 8. Turbidity test for milk.
- 9. Test for pasteurization of milk.
- 10. Coliform test for milk.
- 11. Study of alcoholic fermentation- alcohol from different substractsestimation of percentage of alcohol, total acidity and volatile acidity.
- 12. Study of food-spoilage microorganisms in fresh, canned, fermented food and meat.
- 13. Production and analysis of SCP: Spirulina, yeast, Chlorella, mushroom.
- 14. Bioassay methods for vitamins and amino acids.
- 15. Production of yoghurt using specific starter cultures.
- 16. Demonstration of fermenters.
- 17. Production of pectinase by *Aspergillus niger* using wheat bran coffee pulp using small scale fermenter and its assay.
- 18. Production of λ amylase using *Aspergillus oryzai/ Bacillus lichenforis* using bran in small scale solid state fermentation and its assay.
- 19. Production of microbial polysaccharides and yield estimation.

(B) Principles of Recombinant DNA Technology

- 1. Isolation and purification of plasmid from *Agrobacterium tumefaciens* LBA4404.
- 2. To perform restriction digestion of plant genomic DNA and its visualization.
- 3. To perform ligation of insert into the plasmid vector for construction of recombinant plasmid.
- 4. To check the presence of insert in the recombinant plasmid.
- 5. Preparation of competent cells of *E. coli* (strain DH5 α) using CaCl₂ treatment.
- 6. Transformation of *E. coli* (strain DH5 α) by mobilization of plasmid pBsKs into competent *E. coli* cells.
- 7. Demonstration of presence of GUS gene in plasmid of *Agrobacterium tumefaciens* strain LBA4404.
- 8. Preparation of protein samples for profiling on polyacrylamide gel.
- 9. Method of gel casting and sample loading for protein profiling through SDS-PAGE.
- 10. Running of gel, staining, destaining and analysis of protein profiles using standard protein markers.